
Health o meter®

Your Measure for Good Health™

AUTOMATIC WRIST
**Blood Pressure
Monitor**

Model 7632

IMPORTANT NOTES

Please read this entire manual carefully before using this device.

Healthometer®

Automatic Wrist Blood Pressure Monitor

IMPORTANT SAFETY INSTRUCTIONS

DO NOT DESTROY

Features:

- **One-Touch Easy Operation**
- **Large, Easy-To-Read Digital Display**
- **Two-User Memory Feature**
- **50 Automatic Memory Readings for Each User**
- **MicroStep™ Technology Precisely Controls Deflation for Greater Accuracy**
- **Date and Time Indication**
- **Auto Shut-Off**
- **Convenient Storage Case**
- **Compact and Portable**

In accordance with International Standards:

CE **ISO**
0366 9001

Table of Contents

| | |
|---|-------|
| Introduction | 4 |
| Quick Reference Guide | 5 |
| Standards of Quality | 6 |
| About the Wrist Monitor | 6 |
| About Blood Pressure | 6-7 |
| Blood Pressure Standard | 7 |
| Blood Pressure Fluctuation | 7 |
| Carrying Case | 8 |
| Product Component Identification | 8 |
| Display Explanations | 9 |
| Installing Batteries | 9-10 |
| Setting the Time | 10 |
| HOW TO USE | 11-15 |
| • Helpful Tips | 11 |
| • Applying the Cuff | 11 |
| • Proper Posture During Measurement | 12 |
| • Taking a Measurement | 13-14 |
| • Recalling Memory | 14-15 |
| • Clearing Values from Memory | 15 |
| Troubleshooting | 16 |
| Precautionary Notes | 17 |
| Specifications | 18 |
| Patents | 18 |
| Warranty Information | 19 |
| Blood Pressure Log | 20-21 |

Introduction

Congratulations on your purchase of the new Health o meter® Wrist Blood Pressure Monitor Model 7632. This fully-automatic, wrist-mounted device is both easy to use and ideally suited for daily measurements. The large display includes systolic, diastolic and pulse measurements, all clearly shown at the completion of each reading. With its enhanced, comprehensive features, this Health o meter® monitor will be the only blood pressure monitor you will need for your family.


With this blood pressure monitor, you can store up to 50 measurements into 2 Memory Zones, one for each user. This device is ideal for users who want to monitor and track their blood pressure on a regular basis. The Health o meter® blood pressure monitor is compact and portable, making it ideal for both home and travel usage.

Blood pressure measurements determined with the Health o meter® blood pressure monitor are equivalent to those using the cuff/stethoscope auscultation method and are within the limits prescribed by the American National Standards Institute (ANSI). This monitor is to be used by an adult and is intended for in-home use.

Note: Please read this manual carefully before use. For specific information about your own blood pressure, contact your physician. Please be sure to keep this manual for future reference.

Quick Reference

A. Installing Batteries


Battery Type: 2 AAA 1.5 Volt
(See Page 9)

B. Applying the Cuff


Apply cuff to **LEFT** wrist with
PALM FACING UP.
(See Page 11)

C. Correct Posture


Place the wrist at heart level.
(See Page 12) Relax a few minutes
before starting the Blood Pressure
Monitor. (See Helpful Tips, Page 11)

D. Press ON/OFF/START Button


Remain still during measurement.
This unit will automatically inflate
then deflate. (See Page 13)

E. Automatic Memory

Your new Health o meter® monitor
stores a maximum of 50 measured
values for up to 2 users. The symbol
“mem” will blink when measurement
is completed and stored in memory.
(See Page 14)

F. Recall Memory

Use  **User-Selection** button to
select a Memory Zone, then press “**M**”
to recall prior measured values from
that Memory zone. (See Page 14)

Standards of Quality

This Blood Pressure Monitor conforms with the American National Standards Institute/American Association of Medical Instrumentation Standards for Non-Invasive Blood Pressure Monitors. In addition, the device complies with the European regulations and bears the CE mark "CE 0366". The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC dated 14 June 1993 concerning medical devices, as well as the EMC directive 89/336/EEC:

| | |
|-----------------------|---|
| ANSI/AAMI SB10 | Electronic or automated Sphygmomanometers |
| EN 1060-1 (Dec. 1995) | Non-invasive blood pressure measuring equipment General requirements |
| EN 1060-3 (Feb. 1997) | Non-invasive blood pressure measuring equipment Supplementary requirements for electromechanical blood pressure measuring systems |

About the Wrist Blood Pressure Monitor

Your Health o meter[®] blood pressure monitor uses the oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure, which is equivalent to the air pressure. The measurement of your blood pressure is based on this baseline pressure. Once the cuff inflates to block the passage of blood in the artery, the deflation process begins. During the deflation of the wrist cuff, the unit is detecting the pressure oscillations generated between beats of the pulse. Any muscle movement during this period of time will cause measurement error. After detecting the amplitude and the slope of the pressure oscillations during the deflation process, your Health o meter[®] monitor will determine the systolic and diastolic pressures for you. Your pulse rate is detected and displayed at the same time.

ABOUT BLOOD PRESSURE

What is Blood Pressure?


Blood pressure is the pressure exerted on the inside of the artery wall while blood flows through the artery. The pressure measured when the heart contracts and sends blood out of the heart is the "systolic" reading (higher number). The pressure measured when the heart dilates with blood flowing back into the heart is called the "diastolic" (lowest number) blood pressure.

Why Measure Your Blood Pressure?

Various health problems associated to high blood pressure are quite common with people in today's society. There is a significant correlation of high blood pressure with cardiovascular disease. This correlation has made measuring blood pressure a necessity for those individuals identified with high blood pressure, placing them into a high risk category pertaining to cardiovascular disease. Also, as a person ages, one's blood pressure can change over time, thus it makes good "health sense" to monitor your blood pressure on a regular basis.

Blood Pressure Standard

The World Health Organization (WHO) has developed a blood pressure standard, according to which areas of low and high-risk blood pressure are identified. This standard, however, is a general guideline as individual's blood pressure varies among different people and different age groups, etc.


It is important that you consult with your physician regularly. Your physician can tell you your optimal blood pressure range.

Blood Pressure Fluctuation


Blood pressure fluctuates all the time!

Your blood pressure will change over the course of the month, week and even throughout the day. It is also influenced by seasonal and temperature changes.


Carrying Case

Please store your Health o meter® Blood Pressure Monitor in the carrying case provided as shown below:

Product Component Identification

A. LCD Display

B. Memory Recall Button

C. Wrist Cuff

D. Date/Time Set Button

E. Power Button (for ON/OFF)

F. Battery Cover


G. User Selection Button

H. 2 AAA Size 1.5v Alkaline Batteries
(Included)

Display Explanations

Display:

Memory
Zones

Date/Time Indicator

Systolic Pressure

Diastolic Pressure

Pulse Rate

No. of Memory Sequence

Icons:

Memory Zone 1 : Appears when the monitor is operated in memory zone 1

Memory Zone 2 : Appears when the monitor is operated in memory zone 2

Weak Battery Mark : Appears when batteries should be replaced

Pulse Mark : Shows the pulse rate per minute


Memory Mark : Shows the number of stored measurements

Error Mark : Occurs when an error was made during measurement

Installing Batteries

1. Place your thumb on the side of the battery cover. Pull the battery cover in the direction of the cuff to open.
2. Install lower battery first in the bottom position with the battery ribbon underneath. **Battery Type: 2 Alkaline (AAA) size.**
3. Insert the other battery and pull up the battery ribbon over it.
4. Replace the cover by placing the end tabs of the battery cover in the slots on the top of the battery compartment. Click in the other end to secure the battery cover.

Replace the batteries if:

1. The weak battery mark appears in the display.
2. Nothing appears in the display when the power is switched on.

As the supplied batteries are for test only, they may discharge earlier than batteries you buy in stores. Replace the batteries in pairs.

Remove batteries when unit is not in use for extended periods of time.

Note:  Discard used batteries as directed by your local, environmental codes.

Setting the Time

1. To set/adjust the date and time in the monitor, press the “Set” Button  located on the top surface of the monitor. The display will show a blinking number showing the month.
2. Change the month by pressing the  button. Each press of this button will increase the number by one in a cycling manner. Press the **Set Button**  again to confirm the entry, and the screen will show a blinking number representing the number of the day in the month.
3. Change the date, the hour and the minutes as described in Step 2 above, using the  button to change and the  button to confirm the entries.
4. A “0” will reappear as the Blood Pressure Monitor is ready for measurement again.

Press the Set Button

Month Position Will Blink

Press To Select Month

Press To Confirm Selection

How To Use

Helpful Tips

Here are a few helpful tips to help you obtain more accurate readings:

- Do not measure your blood pressure immediately after consuming a large meal. To obtain more accurate readings, please wait one hour before measuring.
- Do not smoke or drink alcohol before measuring your blood pressure.
- You should not be physically tired or exhausted while taking a measurement.
- It is important that you relax during measurement. Try to take a 15-minute rest before taking a reading.
- Do not take measurements if you are under stress or under tension.
- Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
- If the monitor is stored at very low temperature (near freezing), have it placed at a warm location for at least one hour before using it.
- **Wait about 5 minutes before taking a SECOND blood pressure measurement.**

Applying the Cuff

1. Remove all watches, jewelry, etc., prior to attaching the wrist monitor. Clothing sleeves should be loosely rolled up and the cuff should be wrapped on bare skin for correct measurements.
2. Apply cuff to left wrist with palm facing up as shown in Figure A.
3. Make sure the edge of the cuff is about 3/8 inch from the palm as shown in Figure B.


FIGURE A


FIGURE B

FIGURE C

4. In order to insure accurate measurements, fasten the fabric fastener strap securely around your wrist so there is no extra space between the cuff and the wrist as shown in Figure C.

Note: If the cuff is not wrapped tightly enough, the measurement values may be false. However, do not apply the cuff too tightly as it may constrict the smooth flow of blood to the hand.

5. If your physician has diagnosed you with poor circulation on your left arm, carefully place the cuff around your right wrist as shown in Figure D.

FIGURE D

Proper Posture During Measurement

1. Place your elbow on a table so that the cuff is at the same level as your heart as shown in Figure A.

Note: Your heart is located slightly below your armpit and a bit to the left of the middle of your chest.

Relax your entire body, especially the area between your elbow and fingers.

FIGURE A


FIGURE B

2. If the cuff is not at the same level as your heart or if you can not keep your arm completely still throughout the reading, use a soft object such as a folded towel to support your arm as shown in Figure B. Do not allow hard objects to come in contact with the wrist cuff.
3. Turn your palm upwards.
4. Sit upright in a chair, and take 5-6 deep, relaxing breaths.
5. Avoid leaning back while the measurement is being taken as shown in Figure C.




FIGURE C

Taking a Measurement


1. Press the “**User Selection**” button to select Memory Zone 1 or Memory Zone 2. The measurement that you take will automatically be stored in the Zone selected.

2. Place the cuff on the wrist. Press the “**ON/OFF/START**” button . All digits will light up, checking the display functions. The checking procedure will be completed after about 2 seconds.


3. After all symbols appear, the display will show a blinking “0”. The monitor is “**Ready to Measure**” and will **automatically** inflate the cuff to start measurement.

This monitor will re-inflate automatically to **approximately 220mmHg** if the system detects that your body needs a higher pressure to measure your blood pressure.

- When the measurement is completed, systolic, diastolic, and pulse will be shown simultaneously and be saved automatically into the memory. The measurement is then automatically stored into the pre-designated memory zone. **The last 50 measurement readings will be saved into each memory.**


Note:


- This monitor automatically switches off approximately **1 minute** after the last button operation.
- To interrupt the measurement, simply press any key. The cuff will deflate immediately after a key is pressed.
- During the measurement, do not talk or move your arm or hand muscles.

Recalling Values from Memory


- The monitor has two memory zones (1 and 2). **Each zone can store up to 50 measurements.**

- To recall blood pressure readings from memory, use the **“User-Selection”** button to select the appropriate Memory Zone (1 or 2) from which to recall the readings. Once the proper zone is selected, press the **“M”** button. The memory values will be shown on the LCD display.


- Every new press of the “M” button recalls one set of prior stored values. The most recent memory readings will be recalled first.

Clearing Values from Memory

- Press the “User-Selection”  ↔  button to select Memory Zone 1 or Memory Zone 2.
- Continue to press the “M” button for approximately 5 seconds, then the data in the selected memory zone will be erased automatically.

Press the “M” button for 5 seconds to clear all memory

Note: The data in both Memory Zones (Zone 1 and Zone 2) will be erased if one or both batteries are removed from the battery compartment. *(Please remember to record your previous readings [from both zones] into your Blood Pressure Log prior to replacing the batteries.)* You can also erase the memories stored in both zones by pressing the memory recall key for more than 5 seconds.

Troubleshooting

If a problem should arise during use, please check the following points.

| Symptom | Check Point | Correction |
|---|---|---|
| No display when the POWER button is pressed | Have the batteries run down? | Replace them with 2 new "AAA" alkaline batteries |
| | Have the batteries been properly positioned in the compartment? | Re-insert the batteries in the correct position |
| EE mark shown on the display or the blood pressure value is displayed excessively low (or high) | Is the cuff placed correctly? | Wrap the cuff properly so that it is positioned correctly |
| | Did you talk or move during the measurement? | Measure it again quietly |
| | Did you shake the wrist with the cuff on it? | |

Note: If the unit still does not work, please see the Warranty Information in this manual.
Under no circumstance should you disassemble the unit and attempt repairs yourself.

Precautionary Notes:

1. Your Health o meter[®] monitor contains high-precision components. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping or sharp impacts to the main unit, and protect it from dust.
2. Clean the blood pressure monitor body and the cuff carefully with a slightly dampened soft cloth. Do not press on the unit too hard. Do not wash the cuff or use chemical cleaners on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
3. Leaky batteries can damage the unit. Remove the batteries when the unit is not used for an extended period of time.
4. Do not press the ON/OFF/START key  when the cuff is not placed around the wrist.
5. To avoid hazardous situations, the unit should not be operated by children.
6. Do not disassemble the main unit or cuff.
7. If the unit is stored near a freezing temperature, allow it to acclimate at room temperature before use.
8. Your Health o meter[®] Blood Pressure Monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to make any adjustments to it. If the monitor is not working properly please contact Sunbeam Products, Inc. as outlined under the Warranty information in this manual.
9. For users diagnosed with common arrhythmia (*atrial or ventricular premature beats or atrial fibrillation*), diabetes, poor circulation of blood, kidney problems, or for users who have suffered from stroke, or for unconscious users, this device may not be suitable for use. Please consult your physician prior to use.

Specifications:

| | |
|-----------------------|--|
| Power Source | : 3V DC Two (AA) Alkaline Batteries |
| Measurement Method | : Oscillometric |
| Measurement Range | : Pressure: 20~280mmHg : Pulse: 40~180 beats/ minute |
| Accuracy | : Pressure: \pm 3mmHg : Pulse: \pm 5% of reading |
| Pressure Sensor | : Semiconductor |
| Inflation | : Pump Driven Inflation |
| Deflation | : Automatic Pressure Release Valve |
| Memory capacity | : 50 memories for each of 2 zones |
| Auto-shut-off | : 1 minute after last button operation |
| Operation Environment | : Temperature: 5°C ~ 45°C (41°F~101.3°F); 40 ~ 85% RH max. |
| Storage Environment | : -10°C ~ 60°C (14°F~140°F) 10 ~ 95% RH max. |
| Weight | : 130 g (without batteries) |
| Wrist Circumference | : 13.5 ~ 22 cm (5.3 inches ~ 8.7 inches) |
| Limited patient | : Over 18 years old |
| Dimensions | : 79 (L) x 72 (W) x 71 (H) mm [3.1 in (L) x 2.835 in (W) x 2.795 in (H)] |
| K.B. | : Short time operation is 2 minutes : Classification of type BF |

**To ensure accurate measurements, it is recommended that calibration be completed every two years.*

**Specifications are subject to change without notice.*

Patents:

The Health o meter[®] Model 7632 is patented in countries named below:

| | |
|------------------|------------------------------------|
| Germany : | # 295-11-738.9 # 297-21-574.4 |
| Japan : | # 3020497 |
| Taiwan : | # 140134 |
| Mainland China : | # ZL 95217421.9 # ZL 97249758.7 |

The Health o meter[®] Model 7632 is patent pending in other countries.

Warranty Information

Your Health o meter® Blood Pressure Monitor is covered under the limited Warranty offered by Sunbeam Products, Inc. ("Sunbeam"). Sunbeam warrants that for a period of five years from the date of purchase, this product will be free from defects in material and workmanship. Sunbeam, at its option, will repair or replace this product or any component of the product found to be defective during the warranty period. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. This is your exclusive warranty. This warranty is valid for the original retail purchaser from the date of initial retail purchase and is not transferable. Keep the original sales receipt. Proof of purchase is required to obtain warranty performance. Sunbeam dealers, service centers or retail stores selling Sunbeam Products do not have the right to alter, modify or in any way change the terms and conditions of this warranty.

This warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, disassembly, repair or alteration by anyone other than Sunbeam or an authorized service center. Further, the warranty does not cover Acts of God or Nature, such as fire, flood, hurricane and/or tornado.

Sunbeam shall not be liable for any incidental or consequential damages caused by the breach of any express or implied warranty. Except to the extent prohibited by applicable law, any implied warranty of merchantability or fitness for a particular purpose is limited in duration to the duration of the above warranty. Some states, provinces or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusion may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or province to province.

For more information, or to obtain warranty service, contact Sunbeam Products, Inc. at 800 672-5625 or visit the Health o meter® website at www.healthometer.com

©2002 Sunbeam Products, Inc. All rights reserved.

Health o meter® is a registered trademark of Sunbeam Products, Inc.

and MicroStep™ is a trademark of Sunbeam Products, Inc.

Distributed by Sunbeam Products, Inc., Boca Raton, FL 33431.

Manual printed in China.

BLOOD PRESSURE LOG

Name: _____

My Target Blood Pressure is: _____

I am to call my healthcare practitioner:

if my blood pressure goes above _____ or falls below _____.

if I have the following symptoms: _____

| <u>DATE</u> | <u>TIME</u> | <u>BLOOD PRESSURE</u> | <u>COMMENTS</u> |
|-------------|-------------|-----------------------|-----------------|
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |
| _____ | _____ | ____/____ | _____ |

BLOOD PRESSURE LOG

Name: _____

My Target Blood Pressure is: _____


I am to call my healthcare practitioner:

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©2002 Sunbeam Products, Inc. Todos los derechos reservados
Health o meter® es una marca registrada de Sunbeam Products, Inc.
y MicroStep™ es una marca comercial de Sunbeam Products, Inc.
Distribuida por Sunbeam Products, Inc., Boca Raton, FL 33431.